

ZH3

get back, faster



ZAMAR[®]

FAST RECOVERY THERAPY UNIT

COLD-HOT, AIR COMPRESSION



NO ICE / NO WATER

Therapeutic cooling technology without ice or water inside.

CLOSED & CLEAN SYSTEM

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off liquid.

INNOVATIVE COOLING

Advance technology based on real-time temperature parameters.

STAY FOCUSED ON TEMPERATURE

Set exact temperature you want to reach.

INTERMITTENT COMPRESSION

Improves benefits of cold therapy and healing with 3 levels of AIR compression.

NEW SYSTEM OF COMPRESSION MANAGEMENT

Possibility to select different AIR compression setup.

INSPIRED BY YOUR SHAPE

Complete range of Anatomic wraps designed on your body.

FLEXIBILITY & POWER





FLEXIBILITY & POWER



ZH3

WHEN EVOLUTION NEVER ENDS

PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

FRIENDLY USE

All treatment parameters manageable in one screen.

REPEAT YOUR FAVORITE

Possibility to store up to 12 personal treatments for frequent use

TWO STAGES OF TEMPERATURE

Set your custom temperature and treatment time





Designed and Developed by HTFenix® GmbH

ZAMAR MEDICAL d.o.o.
 Sv. Martin, 6 - 52450 - Vrsar - (HR)
 Tel: +385 (0)52 496 111 - VAT: 96642432252
 info@zamar.care

Features

APPLICATION	COLD - HOT - Compression
TEMPERATURE RANGE	+0°C +40°C (32°F 104°F)
SCREEN	LCD 7" color TouchScreen
WRAPS	Glycol Precharged
Glycol	NON TOX freeze
Air Compression levels	mmHg: 40-60-80
Medical class for Applied Parts	BF
Number of Wraps	up to 2
Storable therapies	1 + 12
DIMENSIONS	cm. 27 x 31 (37 x 37) x 86h inc. 11 x 12 (15 x 15) x 34h
WEIGHT	25 kg. - (55 lb.)

- LCD 7", TFT color touch screen
- Easy Software, Intuitive Menu icons
- Possibility to select two different temperature and treatment time
- Possibility to manage different levels and durations of the AIR compression
- 12 storable programs that can be customized by user
- Totally automated cycle management
- Fast connection of wraps and pipes with Dripless System
- Anatomic Wraps pre charged with glycol
- NON-TOX cooling glycol

www.zamar.care

