

ZT CLINIC

Med 2.0

get back, faster



FAST RECOVERY THERAPY UNIT

COLD - HOT , AIR COMPRESSION



NO ICE / NO WATER

Therapeutic cooling technology without ice or water inside.

CLOSED & CLEAN SYSTEM

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off the liquid.

STAY FOCUSED ON TEMPERATURE

Set exact temperature you want to reach.

WIDE RANGE OF SETTINGS

Set temperature you want from -5°C (23°F) to 45°C (113°F).

INTERMITTENT COMPRESSION

Improves benefits of cold therapy and healing with 3 levels of AIR compression.

NEW SYSTEM OF COMPRESSION MANAGEMENT

Possibility to select different AIR compression setup.

INSPIRED BY YOUR SHAPE

Complete range of Anatomic wraps designed on your body.

POWER multi patients

ZT CLINIC




* configuration with double THOSE (up to 4 Wraps)



POWER & multi patients

ZT CLINIC



WHEN THE POWER IS NOT ENOUGH

PLUG 'N' PLAY RECOVERY

No set up time and preparation needed.

REPEAT YOUR FAVORITE

Possibility to store up to 12 personal treatments for frequent use.

TWO STAGES OF TEMPERATURE

Ability to set 2 different temperatures and time phases.

MORE POWER AND FLEXIBILITY FOR THE NEW DEVICE

4 wraps can be used in the same time for a better result in many treatments.



ZTClinic MG675

ZT CLINIC

Med 2.0



Designed and Delveloped by HTFenix® GmbH - Austria

ZAMAR MEDICAL d.o.o.

Sv. Martin, 6 - 52450 - Vrsar - (HR)

Tel: +385 (0)52 496 111 - VAT: 96642432252

info@zamar.care

Features

APPLICATION	COLD & HOT
TEMPERATURE RANGE	-5°C + 45°C (23°F 113°F)
SCREEN	LCD 10,4" color TouchScreen
WRAPS	Glycol Precharged
Glycol	NON TOX freeze
Air Compression levels	mmHg: 40 - 60 - 80
Working time	up to 300 h.
Medical class for Applied Parts	BF
Number of Wraps	up to 4
Storable therapies	1 + 12
DIMENSIONS	cm. 37 x 55 x 93h inch 15 x 22 x 37h
WEIGHT	43 kg. - (95 lb.)

- LCD 10,4", TFT color touch screen
- Easy software, Intuitive Menu icons
- Possibility to select two different temperature and treatment time
- Possibility to manage different levels and durations of the AIR compression
- 12 storable programs that can be customized by user
- Totally automated cycle management
- Fast connection of wraps and pipes with Dripless System
- Anatomic Wraps pre charge with glycol
- NON-TOX cooling glycol

www.zamar.care

