

# Z-ONE®

## get back, faster



# FAST RECOVERY

## COLD & COMPRESSION THERAPY UNIT



### NO ICE / NO WATER

Therapeutic cooling technology without ice or water inside.

### CLOSED & CLEAN SYSTEM

ZAMAR DRIPLESS SYSTEM preloaded with glycol doesn't need to refill or drain off liquid.

### INNOVATIVE COOLING

Advance technology based on real-time temperature parameters.

### STAY FOCUSED ON TEMPERATURE

Set exact temperature you want to reach.

### INTERMITTENT COMPRESSION

Improves benefits of cold therapy and healing with 3 levels of AIR compression.

### NEW SYSTEM OF COMPRESSION MANAGEMENT

Possibility to select different AIR compression setup.

### INSPIRED BY YOUR SHAPE

Complete range of Anatomic wraps designed on your body.

# LIGHTER & STRONGER





# LIGHTER & STRONGER

**Z-ONE**

**small and powerful than ever**

## **PLUG 'N' PLAY RECOVERY**

No set up time and preparation needed.

## **FRIENDLY USE**

All treatment parameters manageable in one screen.

## **ALWAYS ON THE MOVE**

Light and handy can easily be with you anytime, anywhere.

## **REPEAT YOUR FAVORITE**

Possibility to store your personal treatments for frequent use



Z-ONE®



Designed and Developed by HTFenix® GmbH

**ZAMAR MEDICAL d.o.o.**

Sv. Martin, 6 - 52450 - Vrsar - (HR)  
Tel: +385 (0)52 496 111 - VAT: 96642432252  
info@zamar.care

Features

APPLICATION	COLD - Compression
TEMPERATURE RANGE	+5°C +15°C (41°F 59°F)
SCREEN	LCD 5" color TouchScreen
WRAPS	Glycol Precharged
Glycol	NON TOX freeze
Air Compression levels	mmHg: 40 - 60 - 80
Medical class for Applied Parts	BF
Number of Wraps	up to 2
Storable therapies	1 + 2
DIMENSIONS	cm. 29 x 26 x 30h inc. 11 x 10 x 12h
WEIGHT	9 kg. - (20 lb.)

- LCD 5", TFT color touch screen
- Easy Software, Intuitive Menu icons
- Possibility to select temperature and treatment time
- Possibility to manage different levels and Duration of the AIR compression
- 2 storable programs that can be customized by user
- Totally automated cycle management
- Fast connection of wraps and pipes with Dripless System
- Anatomic Wraps pre charge with glycol
- NON-TOX cooling glycol

www.zamar.care

